

Culligan



Aqua-Cleer®

RO² MD BiO E

DOUBLE PASS REVERSE OSMOSIS

DOUBLE PASS REVERSE OSMOSIS

The Culligan **RO² MD BiO E** is Culligan's most technologically advanced water treatment system for haemodialysis utilizing the latest in bi-osmosis water treatment technology.

The **RO**² **bi-osmosis** system consists of two separate reverse osmosis systems operating in series to produce water of the highest quality and purity.

By operating in series the system effectively treats the water twice. In the first pass, 90-99% of the salts and other contaminants in the water are removed, and in the second pass, again an additional 90-99% of the residual contaminants are removed.

Most importantly this provides a guaranteed barrier to microbial contamination.


An added advantage to the **Culligan bi-osmosis system** is that in the event of an emergency, the system can operate on a single reverse osmosis system and still produce dialysis quality water.

The Culligan RO² will automatically move to this single pass configuration in the unlikely event of a problem, so there is no need for manual intervention.

- CLASS IIb
- FLOWS UP TO 4400 L/H
- ADVANCED SAFETY SYSTEMS
- FULLY ELECTRONIC
- NOTIFIED BODY NO. CE 0434
- 4 MODELS AVAILABLE

Aqua-Cleer® RO² MD BiO E





AVAILABLE MODELS FOR AQUA-CLEER® RO2 MD BIO E

Models	Flow I/h*(GPM)	Installed Power kw (HP)	Dimensions width. x depth. x height. mm (inches)
RO ² MD BiO 1 E	1000 (4.4)	2×3 (2×4 HP)	1800×800×1900 (70 × 31 × 75 inches)
RO ² MD BiO 2 E	1425 (6.27)	2×3 (2×4 HP)	1800×800×1900 (70 × 31 × 75 inches)
RO ² MD BiO 3 E	2500 (11)	2×4 (2×5.4 HP)	1800×800×1900 (70 × 31 × 75 inches)
RO ² MD BiO 4 E	3800 (16.73)	2x4 (2x5.4 HP)	1800x800x1900 (70 x 31 x 75 inches)

^{*} Note : Flow rates shown relate to new membranes operating on properly pretreated feed water of 800 mg/l TDS, 20 °C

All of the Aqua-Cleer® RO2 MD BiO E

models are PLC controlled

Simple and intuitive "touch screen" interface with advanced controls.



Dynamic views

Inlet to Reverse Osmosis conductivity

Conductivity of water produced by 1st stage RO

Conductivity of water produced by 2nd stage RO

Flow rate of RO water produced by 1st stage RO

Flow rate of RO water produced by 2nd stage RO

Flow rate of reject water produced by 1st stage RO

Flow rate of reject water produced by 2nd stage RO

Flow rate of water returning from dialysis loop

Flow rate and consumption rate of artificial kidneys

RO 1st stage operating pressure

RO 2nd stage operating pressure

Inlet to Reverse Osmosis temperature

Distribution loop temperature

REDOX meter reading

Graphic views

Graphical trending of conductivity

Graphical trending of system flow rates

Graphical trending of system temperature

Historical archives of operating data

Immediate display of error messages

Culligan's bi-osmosis water treatment equipment is the result of 30 years of experience in the haemodialysis sector.

The new Bi-Osmosis system is Class IIb Certified and is medical device certified according to

UNI EN ISO 13485, and CE 0434.



WATER QUALITY AND DIALYSATE GUIDELINES & CULLIGAN CERTIFICATES

Chemical contaminants

Category		Parameter	Measure unit	Monitoring place and maximum acceptable level		
				Raw water (yearly frequency)	Dilution water (every 6 months)	
				Max level	Max level	
1) Natural source	Inorganic:	Calcium	mg/L		2	
		Chloride	mg/L	250	50	
		Hydrogen ions	Unità pH	6.5 ÷ 9.5	4.4 ÷ 7.5	
		Fluoride	mg/L	1.5	0.2	
		Magnesium	mg/L		2	
		Potassium	mg/L		2	
		Sodium	mg/L	200	50	
		Sulphate	mg/L	250	50	
	Inorganic:	Mercury	mg/L	0.001	0.001	
0) 4 1		Lead	mg/L	0.01 (0.025 up to 25/12/2013)	0.1	
2) Anthropogenic		Nitrates (as NO ₃)	mg/L	50	2	
		Ammonium	mg/L	0.5	0.2	
	Organic:	Total organo-alogenated compounds	μg/L	30	30	
3) Products used for	Inorganic:	Aluminium	mg/L	0.2	0.01	
potabilisation	morganit.	Total Chlorine	mg/L	0.2	0.1	
4) Byproducts for potabilisation	Inorganic:	Zinc	mg/L		0.1	
treatment	Organic:	Total THM	μg/L	30	30	

Note: these data refer to Italian Pharmacopoeia

Microbiological controls

	Raw water		Treated water		Regular dialysate		Ultrapure dialysate	
Test	Reference value	Frequency	Reference value	Frequency	Reference value	Frequency	Reference value	Frequency
UFC/mL Bacteria at 22 °C	< 100	every 6 months	< 100	monthly	< 100		< 0.1	every 2 months in each display where online treatments are effected
Molds and Yeasts /mL	-	-	< 10	every 6 months	< 10	every 4 months in each display	0	
Endotoxins UI/mL	-	-	< 0.25	monthly	< 0.25		< 0.03	





PLACE YOUR WATER TREATMENT NEEDS IN THE HANDS OF A GLOBAL LEADER



Contact your local Culligan representative for more details

Culligan reserves the right to change any technical or design specifications for the models shown in this brochure.

Place your commercial and industrial water treatment needs in the hands of a global leader.

For over 80 years, Culligan has made better water. Our global network, comprised of 900+ dealers and international licensees in over 90 countries, is dedicated to addressing your water-related problems. As a worldwide leader in water treatment, our sales representatives and service technicians are familiar with the local water conditions in your area. Being global and local position us to deliver customized solutions to commercial and industrial water issues that affect your business and your bottom line.